



## **Copley Hospital, Inc. Laboratory Services Directory**

### **INSTRUCTIONS FOR PACKING SAMPLES FOR TRANSPORT**

1. Place samples in a zip-lock biohazard bag. Make sure that samples with different storage requirements are packaged separately.
2. Make sure the zip lock seal is closed.
3. Put the requisitions, corresponding to the samples, in the back pocket of the bag.
4. Store the specimens as indicated, until delivered to the lab.

#### **NOTES**

- a. Urine and histology specimen cups and microbiology specimens should be individually packaged with the requisitions in the outside pocket.
- b. Make sure that specimen containers (such as urines, feces, and sputums) have their covers securely tightened to prevent leakage.